

For Research Use Only

# FcZero-rAb™ CoraLite® Plus 488 Anti-Mouse CD86 Rabbit Recombinant Antibody

Catalog Number: CL488-FcA98025



## Basic Information

Catalog Number:

CL488-FcA98025

Size:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0802

GenBank Accession Number:

NM\_019388

GeneID (NCBI):

12524

UNIPROT ID:

P42082-1

Full Name:

CD86 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230476B8

Recommended Dilutions:

FC: 0.25 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : LPS treated mouse splenocytes,

## Background Information

CD86 (also known as B7-2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and CTLA-4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with CTLA-4 negatively regulates T-cell activation and diminishes the immune response. (PMID: 7513726; 1847722; 11029388; 27591335)

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3




For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> untreated or LPS-treated mouse splenocytes were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse CD86 Rabbit RecAb (CL488-FcA98025, Clone: 230476B8). Cells were not fixed.

1x10<sup>6</sup> LPS-treated mouse splenocytes were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse CD86 Rabbit RecAb (CL488-FcA98025, Clone: 230476B8) (red) or 0.25 ug CoraLite® Plus 488 Rabbit IgG Isotype Control RecAb (CL488-98136, Clone: 240953C9) (blue). 1x10<sup>6</sup> untreated mouse splenocytes were surface stained with 0.25 ug CoraLite® Plus 488 Anti-Mouse CD86 Rabbit RecAb (CL488-FcA98025, Clone: 230476B8) (black, dashed). Cells were not

