

For Research Use Only

FcZero-rAb® CoraLite® Plus 488 Anti-Human IL-2 Rabbit Recombinant Antibody

Catalog Number: **CL488-FcA98053**



Basic Information

Catalog Number:

CL488-FcA98053

Size:

100tests , 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0828

GenBank Accession Number:

BC070338

GeneID (NCBI):

3558

UNIPROT ID:

P60568

Full Name:

interleukin 2

Calculated MW:

153 aa, 18 kDa

Purification Method:

Protein A purification

CloneNo.:

240416E3

Recommended Dilutions:

FC (Intra): 5 ul per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Positive Controls:

FC (Intra) : PMA and Iono treated human PBMCs,

Background Information

Interleukin-2 (IL-2) is a four-helix bundle, type I cytokine that functions as a growth factor for a wide range of leukocytes. In the immune system, IL-2 is essential for immune homeostasis, normal T regulatory cell function, and self-tolerance. It regulates immune cell homeostasis and has been used to treat a range of disorders including cancer and autoimmune disease. IL-2 signals through heterodimerization of the IL-2R β and IL-2R γ receptor subunits (PMID: 25992858, 24907378, 15034008, 17981641).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 untreated or PMA and Iono treated human PBMCs were intracellularly stained with FcZero-rAb® CoraLite® Plus 750 Anti-Human CD3 and 5 μ l CoraLite® Plus 488 Anti-Human IL-2 Rabbit RecAb (CL488-FcA98053, Clone: 240416E3). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).

1×10^6 PMA and Iono treated human PBMCs were intracellularly stained with 5 μ l CoraLite® Plus 488 Anti-Human IL-2 Rabbit RecAb (CL488-FcA98053, Clone: 240416E3) (red) or FcZero-rAb® CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-FcA98136, Clone: 240953C9) (blue). 1×10^6 untreated human PBMCs were intracellularly stained with 5 μ l CoraLite® Plus 488 Anti-Human IL-2 Rabbit RecAb (CL488-FcA98053, Clone: 240416E3)

