

For Research Use Only

FcZero-rAb® CoraLite® Plus 488 Anti-Human TREM-1/CD354 Rabbit Recombinant Antibody

Catalog Number: CL488-FcA98154



Basic Information

Catalog Number:

CL488-FcA98154

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1438

GenBank Accession Number:

NM_018643.5

GeneID (NCBI):

54210

UNIPROT ID:

Q9NP99-1

Full Name:

triggering receptor expressed on myeloid cells 1

Calculated MW:

26kDa

Purification Method:

Protein A purification

CloneNo.:

241279D10

Recommended Dilutions:

FC: 5 ul per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : human peripheral blood leukocytes,

Background Information

Triggering receptor expressed on myeloid cells 1 (TREM-1, also known as CD354) is a member of the immunoglobulin superfamily, which is expressed on neutrophils and monocytes. TREM1 plays a role in the innate immune system. TREM1 amplifies neutrophil and monocyte-mediated inflammatory responses triggered by bacterial and fungal infections by stimulating release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were surface stained with FcZero-rAb® CoraLite® Plus 405 Anti-Human CD16 (3G8), and 5 ul CoraLite® Plus 488 Anti-Human TREM-1/CD354 Rabbit RecAb (CL488-FcA98154, Clone: 241279D10) or FcZero-rAb® CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10⁶ human peripheral blood leukocytes were surface stained with 5 ul CoraLite® Plus 488 Anti-Human TREM-1/CD354 Rabbit RecAb (CL488-FcA98154, Clone: 241279D10) (red) or FcZero-rAb® CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.