

For Research Use Only

FcZero-rAb® CoraLite® Plus 488 Anti-Human CD40L/CD154 Rabbit Recombinant Antibody

Catalog Number: CL488-FcA98355-2



Basic Information

Catalog Number:
CL488-FcA98355-2

Size:
100 tests, 5 µl/test

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
EG3501

GenBank Accession Number:
BC071754

GeneID (NCBI):
959

UNIPROT ID:
P29965

Full Name:
CD40 ligand

Calculated MW:
261 aa, 29 kDa

Purification Method:
Protein A purification

CloneNo.:
242847F7

Recommended Dilutions:
FC: 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:
493 nm / 522 nm

Applications

Tested Applications:
FC

Species Specificity:
human

Positive Controls:

FC : PMA and ionomycin treated human PBMCs,

Background Information

The CD40 ligand (CD40L, TRAP, CD154), a member of the TNF superfamily of ligands, is expressed as either a 33 kDa transmembrane homologue or 18 kDa soluble form (sCD154). CD40L is primarily expressed on activated CD4⁺ T cells and on a small proportion of CD8⁺ T cells and platelets. It binds to CD40 on antigen-presenting cells (APC), which leads to many effects depending on the target cell type. Recent studies have suggested that CD40/CD40L interactions regulate oxidative stress and affect various signaling pathways in both the immunological and the cardiovascular systems. The CD40/CD40L system is also involved in tumorigenesis.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or PMA and ionomycin treated human PBMCs were surface stained with FcZero-rAb® CoraLite® Plus 647 Anti-Human CD4 and 5 ul CoraLite® Plus 488 Anti-Human CD40L/CD154 Rabbit RecAb (CL488-FcA98355-2, Clone: 242847F7). Cells were not fixed.

1x10⁶ PMA and ionomycin treated human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD40L/CD154 Rabbit RecAb (CL488-FcA98355-2, Clone: 242847F7) (red) or FcZero-rAb® CoraLite® Plus 488 Rabbit IgG Isotype Control Recombinant Antibody (CL488-FcA98136, Clone: 240953C9) (blue). 1x10⁶ untreated human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD40L/CD154 Rabbit RecAb