

For Research Use Only

CoraLite®555 Mouse IgG2b Isotype Control (MPC-11)



Catalog Number:CL550-65128

Basic Information

Catalog Number: CL550-65128	GenBank Accession Number: GeneID (NCBI): 16016	Purification Method: Affinity purification
Size: 100ug, 200 µg/mL	Full Name: immunoglobulin heavy chain 3 (serum IgG2b)	CloneNo.: MPC-11
Source: Mouse		Excitation/Emission maxima wavelengths: 554 nm / 570 nm
Isotype: IgG2b		

Applications

Tested Applications:
FC
Species Specificity:

Background Information

The MPC-11 immunoglobulin is useful as an isotype-matched control. The MPC-11 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2b antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data