

For Research Use Only

CoraLite®555-conjugated CRCP Polyclonal antibody

Catalog Number:CL555-14348



Basic Information

Catalog Number:

CL555-14348

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5689

GenBank Accession Number:

BC040107

GeneID (NCBI):

27297

UNIPROT ID:

O75575

Full Name:

CGRP receptor component

Calculated MW:

17 kDa

Observed MW:

20 kDa

Purification Method:

Antigen Affinity Purified

Excitation/Emission maxima

wavelengths:

557 nm / 570 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Raji cells were intracellularly stained with 0.5 ug CoraLite®555 Anti-Human CRCP (CL555-14348) (red), or 0.5 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).