

For Research Use Only

CoraLite®555-conjugated CRCP Polyclonal antibody



Catalog Number:CL555-14348

Basic Information

Catalog Number: CL555-14348	GenBank Accession Number: BC040107	Purification Method: Antigen Affinity Purified
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 27297	Excitation/Emission maxima wavelengths: 555 nm / 570 nm
Source: Rabbit	Full Name: CGRP receptor component	
Isotype: IgG	Calculated MW: 17 kDa	
Immunogen Catalog Number: AG5689	Observed MW: 20 kDa	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Raji cells were intracellularly stained with 0.5 ug CoraLite®555 Anti-Human CRCP (CL555-14348) (red), or 0.5 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).