

For Research Use Only

CoraLite®555 Anti-Mouse CD3 (17A2)



Catalog Number: **CL555-65077**

Basic Information

Catalog Number: CL555-65077	GenBank Accession Number: BC098236	Purification Method: Affinity purification
Size: 100ug, 200 µg/ml	GeneID (NCBI): 12501	CloneNo.: 17A2
Source: Rat	Full Name: CD3 antigen, epsilon polypeptide	Excitation/Emission maxima wavelengths: 555 nm / 570 nm
Isotype: IgG2b, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

CD3 is a multimeric protein associated with the T-cell receptor (TCR) to form a complex involved in antigen recognition and signal transduction (PMID: 15885124). CD3 is composed of CD3γ, δ, ε, and ζ chains (PMID: 1826255). It is expressed by thymocytes in a developmentally regulated manner, T cells, and some NK cells (PMID: 3289580). The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction (PMID: 11985657).

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with CoraLite® Plus 647 Anti-Mouse CD4 (GK1.5) (CL647-65104, Clone: GK1.5) and 0.5 ug CoraLite® 555 Anti-Mouse CD3 (CL555-65077, Clone: 17A2) or 0.5 ug Isotype Control. Cells were not fixed.