

For Research Use Only

CoraLite® Plus 555 Mouse IgG2a Isotype Control (C1.18.4)



Catalog Number: CL555-65208

Basic Information

Catalog Number:
CL555-65208

Size:
100ug, 0.2 mg/ml

Source:
Mouse

Isotype:
IgG2a, kappa

GenBank Accession Number:
GeneID (NCBI):

Full Name:

Purification Method:
Affinity purification

CloneNo.:
C1.18.4

Excitation/Emission maxima
wavelengths:
554 nm / 570 nm

Applications

Tested Applications:
FC

Species Specificity:

Background Information

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 0.2 ug CoraLite® Plus 555 Anti-Mouse IgG2a Isotype Control (CL555-65208, Clone:C1.18.4) (red), or CL555 Anti-Human CD3 (OKT3) (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.