

For Research Use Only

CoraLite® Plus 555 Anti-Mouse CD48 (HM48-1)



Catalog Number: **CL555-65227**

Basic Information

Catalog Number: CL555-65227	GenBank Accession Number: BC060977	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	GeneID (NCBI): 12506	CloneNo.: HM48-1
Source: Armenian Hamster	Full Name: CD48 antigen	Excitation/Emission maxima wavelengths: 554 nm / 570 nm
Isotype: IgG		

Applications

Tested Applications:
FC
Species Specificity:
Mouse

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 555 Anti-Mouse CD48 (CL555-65227, Clone:HM48-1) (red), or 0.5 ug Isotype Control. Cells were not fixed.