

For Research Use Only

CoraLite® Plus 555 Anti-Rat CD43 (W3/13)



Catalog Number: **CL555-65237**

Basic Information

Catalog Number:

CL555-65237

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

24796

UNIPROT ID:

P13838

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

N/A

CloneNo.:

W3/13

Excitation/Emission maxima wavelengths:

554 nm / 570 nm

Applications

Tested Applications:

FC

Species Specificity:

Rat

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 555 Anti-Rat CD43 (CL555-65237, Clone:W3/13), or 0.5 ug CoraLite® 555 Mouse IgG1 Isotype Control (MOPC-21) (CL555-65124, Clone: MOPC-21). Cells were not fixed.