

For Research Use Only

# CoraLite®555-conjugated Fascin Monoclonal antibody



Catalog Number:CL555-66321

## Basic Information

<b>Catalog Number:</b> CL555-66321	<b>GenBank Accession Number:</b> BC000521	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 6624	<b>CloneNo.:</b> 1B5G3
<b>Source:</b> Mouse	<b>Full Name:</b> fascin homolog 1, actin-bundling protein (Strongylocentrotus purpuratus)	<b>Recommended Dilutions:</b> IF 1:100-1:500
<b>Isotype:</b> IgG2a	<b>Calculated MW:</b> 55 kDa	<b>Excitation/Emission maxima wavelengths:</b> 557 nm/570 nm
<b>Immunogen Catalog Number:</b> AG6074	<b>Observed MW:</b> 55 kDa	

## Applications

<b>Tested Applications:</b> FC (Intra), IF	<b>Positive Controls:</b> IF : human tonsillitis tissue,
<b>Species Specificity:</b> human	

## Background Information

Fascin is a highly conserved and ubiquitous actin cross-linking protein that has a major function in cell motility and adhesion. Fascin localizes to a number of highly dynamic cellular structures that require strong mechanical support which include stress fibers and cellular protrusions such as microvilli, microspikes and lamellipodia. Fascin is up-regulated in many human carcinomas. It's a marker for Reed-Sternberg cells in Hodgkin's disease; it is also a highly selective marker for dendritic cells of lymphoid tissues.



## Storage

**Storage:**  
Store at -20°C. Avoid exposure to light.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> HeLa cells were intracellularly stained with 0.4 ug CoraLite®555 Anti-Human Fascin (CL555-66321, Clone:1B5G3) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).

Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®555 Fascin antibody (CL555-66321, Clone: 1B5G3) at dilution of 1:200.