

For Research Use Only

# CoraLite®555 Mouse IgG2b isotype control (11B8C4)



Catalog Number:CL555-66360-3

## Basic Information

Catalog Number:  
CL555-66360-3

GenBank Accession Number:  
GeneID (NCBI):

Purification Method:  
Protein A purification

Size:  
100ug, Concentration: 0.2 mg/mL by  
Nanodrop;

Full Name:

CloneNo.:  
11B8C4

Source:  
Mouse

Excitation/Emission maxima  
wavelengths:  
557 nm / 570 nm

Isotype:  
IgG2b

## Applications

Tested Applications:  
FC

Species Specificity:  
Megathura crenulata

## Background Information

This antibody is used as the isotype control of mouse IgG2b. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

## Storage

Storage:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.09% sodium azide and 0.5% BSA.

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite®555 Anti-Mouse IgG2b isotype control (CL555-66360-3, Clone:11B8C4) (red), or CoraLite®555 Anti-Human CD4 (OKT4) (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.