

For Research Use Only

CoraLite® Plus 555-conjugated NSFL1C Monoclonal antibody

Catalog Number:CL555-68182



Basic Information

Catalog Number: CL555-68182	GenBank Accession Number: BC002801	Purification Method: Protein G purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 55968	CloneNo.: 1F12E12
Source: Mouse	UNIPROT ID: Q9UNZ2	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG1	Full Name: NSFL1 (p97) cofactor (p47)	Excitation/Emission maxima wavelengths: 554 nm / 570 nm
Immunogen Catalog Number: AG8040	Calculated MW: 41 kDa	
	Observed MW: 41 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : A431 cells,
Species Specificity: Human, Mouse, Rat, Pig	


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A431 cells using CoraLite® Plus 555 NSFL1C antibody (CL555-68182, Clone: 1F12E12) at dilution of 1:200.