

For Research Use Only

CoraLite® Plus 555-conjugated ZO-1 Recombinant monoclonal antibody

Catalog Number:CL555-82870-7



Basic Information

Catalog Number: CL555-82870-7	GenBank Accession Number: BC111712	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 7082	CloneNo.: 240150D7
Source: Rabbit	UNIPROT ID: Q07157	Recommended Dilutions: IF/ICC: 1:50-1:500
Isotype: IgG	Full Name: tight junction protein 1 (zona occludens 1)	Excitation/Emission maxima wavelengths: 554 nm / 570 nm
Immunogen Catalog Number: AG33182	Calculated MW: 1748 aa, 195 kDa	
	Observed MW: 230 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : MCF-7 cells,
Species Specificity: human, mouse, canine	

Background Information

Tight junction (or zonula occludens) form the continuous intercellular barrier between epithelial and endothelial cells, which is required to separate tissue spaces and regulate selective movement of solutes across the epithelium and endothelium (PMID: 20066090). ZO-1 (also known as TJP1) is a peripheral membrane phosphoprotein located on the cytoplasmic face and is expressed in tight junctions of both epithelial and endothelial cells (PMID: 3528172). It binds the transmembrane proteins occludin and the claudins linking them to cytoskeletal actin (PMID: 17418867). ZO-1 belongs to a family of multidomain proteins known as the membrane-associated guanylate kinase homologs (MAGUKs). It is a pivotal tight junction protein and may be involved in signaling mechanisms regulating cell proliferation and differentiation (PMID: 22782886).


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite® Plus 555 ZO-1 antibody (CL555-82870-7, Clone: 240150D7) at dilution of 1:200.