

For Research Use Only

CoraLite®594-conjugated BOULE Polyclonal antibody

Catalog Number:CL594-13720



Basic Information

Catalog Number:

CL594-13720

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4662

GenBank Accession Number:

BC033674

GeneID (NCBI):

66037

UNIPROT ID:

Q8N9W6

Full Name:

bol, boule-like (Drosophila)

Calculated MW:

283 aa, 31 kDa

Observed MW:

38-42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human, mouse, rat

Positive Controls:

IF-P : mouse testis tissue,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite® 594 BOULE antibody (CL594-13720) at dilution of 1:100.



Immunofluorescent analysis of (4% PFA) fixed mouse testis tissue using CoraLite® 594 BOULE antibody (CL594-13720) at dilution of 1:100.