

For Research Use Only

CoraLite®594-conjugated RBP7 Polyclonal antibody



Catalog Number:CL594-14541

Basic Information

Catalog Number: CL594-14541	GenBank Accession Number: BC063013	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 116362	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	Full Name: retinol binding protein 7, cellular	Excitation/Emission maxima wavelengths: 593 nm / 614 nm
Isotype: IgG	Calculated MW: 15.5 kDa	
Immunogen Catalog Number: AG6045	Observed MW: 15-18 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : human kidney tissue,
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human kidney tissue using CoraLite® 594 RBP7 antibody (CL594-14541) at dilution of 1:200.