

For Research Use Only

CoraLite®594-conjugated MOBKL3 Polyclonal antibody

Catalog Number:CL594-15886



Basic Information

Catalog Number:

CL594-15886

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8671

GenBank Accession Number:

BC005237

GeneID (NCBI):

25843

UNIPROT ID:

Q9Y3A3

Full Name:

MOB1, Mps One Binder kinase activator-like 3 (yeast)

Calculated MW:

225 aa, 26 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF-P: 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human, mouse, rat

Positive Controls:

IF-P : human colon cancer tissue,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human colon cancer tissue using CoraLite®594 MOBKL3 antibody (CL594-15886) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).