

For Research Use Only

CoraLite® Plus 594 Anti-Human CD44 (F10-44-2)



Catalog Number: **CL594-65063**

Basic Information

Catalog Number: CL594-65063	GenBank Accession Number: BC004372	Purification Method: Affinity purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 960	CloneNo.: F10-44-2
Source: Mouse	Full Name: CD44 molecule (Indian blood group)	Excitation/Emission maxima wavelengths: 594 nm / 615 nm
Isotype: IgG2a, kappa	Calculated MW: 742 aa, 82 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

CD44 is a type I transmembrane glycoprotein expressed on embryonic stem cells and in various levels on other cell types including connective tissues and bone marrow. CD44 expression is also upregulated in subpopulations of cancer cells and is recognized as a molecular marker for cancer stem cells (PMID: 29747682). It is a cell-surface receptor that mediates cell-cell and cell-matrix interactions through its affinity for hyaluronic acid (HA) and possibly also through its affinity for other ligands (PMID: 10694938). Adhesion with HA plays an important role in cell migration, tumor growth and progression. CD44 is also involved in lymphocyte activation, recirculation and homing, and in hematopoiesis.

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 594 Anti-Human CD44 (CL594-65063, Clone:F10-44-2), or CoraLite®594 Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.



1x10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 594 Anti-Human CD44 (CL594-65063, Clone:F10-44-2), or CoraLite®594 Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.