

For Research Use Only

CoraLite®594 Anti-Mouse CD24 (M1/69)



Catalog Number:CL594-65121

Basic Information

| | | |
|--------------------------------|---------------------------------------|---|
| Catalog Number: CL594-65121 | GenBank Accession Number: BC075622 | Purification Method: Affinity purification |
| Size: 100ug, 0.5 mg/ml | GeneID (NCBI): 12484 | CloneNo.: M1/69 |
| Source: Rat | Full Name: CD24a antigen | Excitation/Emission maxima wavelengths: 588 nm / 604 nm |
| Isotype: IgG2b, kappa | | |

Applications

Tested Applications:
FC
Species Specificity:
Mouse

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 0.5 ug CoraLite® 594 Anti-Mouse CD24 (CL594-65121, Clone:M1/69) (red), or 0.5 ug Isotype Control. Cells were not fixed.