

For Research Use Only

CoraLite®594 Rat IgG2a Isotype Control (2A3)



Catalog Number:CL594-65209

Basic Information

Catalog Number:

CL594-65209

Size:

100ug , 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Unknown

CloneNo.:

2A3

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

Background Information

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 BALB/c mouse splenocytes were surface stained with 0.5 μ g CoraLite®594 Anti-Rat IgG2a Isotype Control (CL594-65209, Clone: 2A3) (red), or 0.5 μ g CoraLite®594 Anti-Mouse CD16/32 (93) (CL594-65057, Clone: 93) (green) or unstained (blue). Cells were not fixed.