

For Research Use Only

# CoraLite®594 Anti-Rat CD43 (W3/13)



Catalog Number: **CL594-65237**

---

## Basic Information

<b>Catalog Number:</b> CL594-65237	<b>GenBank Accession Number:</b> GeneID (NCBI): 24796	<b>Purification Method:</b> N/A
<b>Size:</b> 100ug, 0.5 mg/ml	<b>Full Name:</b> sialophorin	<b>CloneNo.:</b> W3/13
<b>Source:</b> Mouse	<b>Calculated MW:</b> 38 kDa	<b>Excitation/Emission maxima wavelengths:</b> 588 nm / 604 nm
<b>Isotype:</b> IgG1, kappa		

---

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Rat

---

## Background Information

---

### Storage

**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.


---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> rat splenocytes cells were surface stained with 0.5 ug CoraLite® 594 Anti-Rat CD43 (CL594-65237, Clone: W3/13), or 0.5 ug Isotype Control. Cells were not fixed.