

For Research Use Only

# CoraLite® Plus 594 Mouse IgG2a Isotype Control (MOPC-173)



Catalog Number: **CL594-65244**

## Basic Information

Catalog Number:

**CL594-65244**

Size:

**100ug , 0.2 mg/ml**

Source:

**Mouse**

Isotype:

**IgG2a, kappa**

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

**N/A**

CloneNo.:

**MOPC-173**

Excitation/Emission maxima  
wavelengths:

**594 nm / 615 nm**

## Applications

Tested Applications:

**FC**

Species Specificity:

## Storage

Storage:

**Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.**

Storage Buffer:

**PBS with 0.09% sodium azide and 0.5% BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite® Plus 594 Anti-Mouse IgG2a Isotype Control (CL594-65244, Clone:MOPC-173) (red), unstained (blue), or CL594 Anti-Human CD3 (Hit3a) (green). Cells were not fixed. Lymphocytes were gated.