

For Research Use Only

CoraLite® Plus 594 Mouse IgG2a Isotype Control (MOPC-173)



Catalog Number: CL594-65244

Basic Information

Catalog Number:
CL594-65244

Size:
100ug, 0.2 mg/ml

Source:
Mouse

Isotype:
IgG2a, kappa

GenBank Accession Number:
GeneID (NCBI):

Full Name:

Purification Method:
N/A

CloneNo.:
MOPC-173

Excitation/Emission maxima
wavelengths:
594 nm / 615 nm

Applications

Tested Applications:
FC

Species Specificity:

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



CL594-65244 Mouse IgG2a (MOPC-173)

1X10⁶ human PBMCs were surface stained with 0.2 ug CoraLite® Plus 594 Anti-Mouse IgG2a Isotype Control (CL594-65244, Clone:MOPC-173) (red), unstained (blue), or CL594 Anti-Human CD3 (Hit3a) (green). Cells were not fixed. Lymphocytes were gated.