

For Research Use Only

CoraLite® Plus 594 Anti-Human CD96 (NK92.39)

Catalog Number: CL594-65250



Basic Information

Catalog Number:

CL594-65250

Size:

100 tests, 5 ul/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC020749

GeneID (NCBI):

10225

UNIPROT ID:

P40200

Full Name:

CD96 molecule

Calculated MW:

402 aa, 46 kDa

Purification Method:

N/A

CloneNo.:

NK92.39

Excitation/Emission maxima wavelengths:

594 nm / 615 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ PHA treated human PBMCs were surface stained with CL488 Anti-Human CD19 (CL488-65197, Clone:4G7) and 5 ul CoraLite® Plus 594 Anti-Human CD96 (CL594-65250, Clone:NK92.39) or CoraLite® 594 Mouse IgG1 Isotype Control (CL594-65124, Clone: MOPC-21). Cells were not fixed. Lymphocytes were gated.