

For Research Use Only

CoraLite®594-conjugated Glutamine Synthetase Monoclonal antibody



Catalog Number:CL594-66323-2

Featured Product

Basic Information

Catalog Number:

CL594-66323-2

Size:

100ul , Concentration: 1000 µg/ml by

Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6309

GenBank Accession Number:

BC011700

GeneID (NCBI):

2752

UNIPROT ID:

P15104

Full Name:

glutamate-ammonia ligase
(glutamine synthetase)

Calculated MW:

374 aa, 42 kDa

Observed MW:

42 kDa

Purification Method:

Protein G purification

CloneNo.:

1B6G4

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

