

For Research Use Only

CoraLite®594-conjugated ANGPTL8 Monoclonal antibody

Catalog Number:CL594-66641



Basic Information

Catalog Number:

CL594-66641

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20708

GenBank Accession Number:

BC146552

GeneID (NCBI):

55908

UNIPROT ID:

Q6UXH0

Full Name:

hepatocellular carcinoma-associated
gene TD26

Calculated MW:

251 aa, 28 kDa

Purification Method:

Protein G purification

CloneNo.:

2H9F9

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima

wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Species Specificity:

Human

Positive Controls:

IF-P : human liver tissue,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human liver tissue using Coralite® 594 ANGPTL8 antibody (CL594-66641, Clone: 2H9F9) at dilution of 1:200.