

For Research Use Only

# CoraLite®594-conjugated INTS3 Monoclonal antibody



Catalog Number:CL594-67060

## Basic Information

<b>Catalog Number:</b> CL594-67060	<b>GenBank Accession Number:</b> BC025254	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 2000 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 65123	<b>CloneNo.:</b> 2G11B7
<b>Source:</b> Mouse	<b>Full Name:</b> integrator complex subunit 3	<b>Recommended Dilutions:</b> IF 1:50-1:500
<b>Isotype:</b> IgG2b	<b>Calculated MW:</b> 118 kDa	<b>Excitation/Emission maxima wavelengths:</b> 593 nm / 614 nm
<b>Immunogen Catalog Number:</b> AG9917	<b>Observed MW:</b> 118 kDa	

## Applications

<b>Tested Applications:</b> IF	<b>Positive Controls:</b> IF : HeLa cells,
<b>Species Specificity:</b> Human, Mouse, Rat	

## Background Information

### Storage


**Storage:**  
Store at -20°C. Avoid exposure to light.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CL594-67060 (INTS3 antibody) at dilution of 1:100.