

For Research Use Only

CoraLite®594-conjugated ARL2BP Monoclonal antibody

Catalog Number:CL594-67316



Basic Information

Catalog Number:

CL594-67316

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29261

GenBank Accession Number:

BC003087

GeneID (NCBI):

23568

UNIPROT ID:

Q9Y2Y0

Full Name:

ADP-ribosylation factor-like 2 binding
protein

Calculated MW:

19 kDa

Observed MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1G5C2

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite®594 Anti-Human ARL2BP (CL594-67316, Clone:1G5C2) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).