

For Research Use Only

CoraLite®594-conjugated TYK2 Monoclonal antibody



Catalog Number:CL594-67411

Basic Information

Catalog Number: CL594-67411	GenBank Accession Number: BC014243	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 7297	CloneNo.: 1E1A10
Source: Mouse	Full Name: tyrosine kinase 2	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2b	Calculated MW: 1187 aa, 134 kDa	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG9442	Observed MW: 134 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : HepG2 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite®594 TYK2 antibody (CL594-67411, Clone: 1E1A10) at dilution of 1:200.