

For Research Use Only

CoraLite®594-conjugated FES Monoclonal antibody

Catalog Number:CL594-67524



Basic Information

Catalog Number:

CL594-67524

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28108

GenBank Accession Number:

BC035357

GeneID (NCBI):

2242

UNIPROT ID:

P07332

Full Name:

feline sarcoma oncogene

Calculated MW:

822 aa, 94 kDa

Observed MW:

93 kDa, 79 kDa, 87 kDa, 85 kDa

Purification Method:

Protein A purification

CloneNo.:

1E1B11

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human

Positive Controls:

IF/ICC : A549 cells,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol)
fixed A549 cells using Coralite®594 FES antibody
(CL594-67524, Clone: 1E1B11) at dilution of 1:200.