

For Research Use Only

CoraLite®594-conjugated HNRNPUL2-Specific Recombinant antibody

Catalog Number:CL594-85739



Basic Information

Catalog Number:

CL594-85739

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001079559

GeneID (NCBI):

221092

UNIPROT ID:

Q1KMD3

Full Name:

heterogeneous nuclear ribonucleoprotein U-like 2

Calculated MW:

85 kDa

Purification Method:

Protein A purification

CloneNo.:

250029G10

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® 594 HNRNPUL2-Specific antibody (CL594-85739, Clone: 250029G10) at dilution of 1:200, CL488-Phalloidin (green).