

For Research Use Only

CoraLite®594 Anti-Human LIGHT/CD258 Rabbit Recombinant Antibody

Catalog Number: CL594-98063



Basic Information

Catalog Number:

CL594-98063

Size:

100tests , 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG32095

GenBank Accession Number:

NM_003807

GeneID (NCBI):

8740

UNIPROT ID:

O43557

Full Name:

tumor necrosis factor (ligand)
superfamily, member 14

Calculated MW:

26KD

Purification Method:

Protein A purification

CloneNo.:

240694D7

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : PMA and ionomycin treated human PBMCs,

Background Information

LIGHT, also known as TNFSF14, HVEM, or CD258, is a type II transmembrane protein that belongs to the TNF superfamily (PMID: 9462508). It is predominantly expressed in lymphoid tissues and produced by activated T cells and immature dendritic cells (DCs) (PMID: 9462508; 10754304). Through its two primary functional receptors HVEM (TNFSF14) and lymphotoxin β receptor (LT β R), LIGHT signaling is involved in development and maintenance of lymphoid tissues, as well as innate and adaptive immune responses (PMID: 24575096; 26720335). In addition to the membrane form, LIGHT can also exist as a soluble form generated by cleavage of the extracellular portion of the membrane form.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or PMA, ionomycin and Brefeldin A treated human PBMCs were surface stained with APC Anti-Human CD8a and 5 ul CoraLite®594 Anti-Human LIGHT/CD258 Rabbit RecAb (CL594-98063, Clone: 240694D7). Cells were not fixed.

1x10⁶ PMA, ionomycin and Brefeldin A treated human PBMCs were surface stained with 5 ul CoraLite®594 Anti-Human LIGHT/CD258 Rabbit RecAb (CL594-98063, Clone: 240694D7) (red) or CoraLite®594 Rabbit IgG Isotype Control RecAb (CL594-98136, Clone: 240953C9) (blue). 1x10⁶ untreated human PBMCs were surface stained with 5 ul CoraLite®594 Anti-Human LIGHT/CD258 Rabbit RecAb (CL594-98063, Clone: 240694D7) (black, dashed). Cells

