

For Research Use Only

CoraLite®594 Anti-Human Granzyme B Rabbit Recombinant Antibody

Catalog Number: CL594-98070



Basic Information

Catalog Number:

CL594-98070

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC030195

GeneID (NCBI):

3002

UNIPROT ID:

P10144

Full Name:

granzyme B (granzyme 2, cytotoxic T-lymphocyte-associated serine esterase 1)

Calculated MW:

247 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

240761G4

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Background Information

GZMB (Granzyme B) is also named as CGL1, CSPB, CTLA1, GRB and belongs to the Granzyme subfamily. This enzyme is necessary for target cell lysis in cell-mediated immune responses. The cytotoxic lymphocyte protease granzyme B (GzmB) can promote apoptosis through direct processing and activation of members of the caspase family. GzmB can also cleave the BH3-only protein, Bcl-2, to promote caspase-independent mitochondrial permeabilization (PMID:17283187). GzmB induces lamin B degradation in isolated nuclei less efficiently than GzmA (PMID:11331782). This full-length protein has 2 glycosylation sites and a signal peptide. Unglycosylated human granzyme B is 26 kDa and high mannose glycosylated is 32 kDa and only 32 kDa or smaller forms of granzyme B are accumulated within nuclei (PMID:8626751). GzmB also forms dimers.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were intracellularly stained with FITC Plus Anti-Human CD8a, and 5 ul CoraLite®594 Anti-Human Granzyme B Rabbit RecAb (CL594-98070, Clone: 240761G4) or CoraLite®594 Rabbit IgG Isotype Control RecAb (CL594-98136, Clone: 240953C9). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).



1x10⁶ human PBMCs were intracellularly stained with 5 ul CoraLite®594 Anti-Human Granzyme B Rabbit RecAb (CL594-98070, Clone: 240761G4) (red), or CoraLite®594 Rabbit IgG Isotype Control RecAb (CL594-98136, Clone: 240953C9) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).