

For Research Use Only

CoraLite®594 Mouse IgG2a, λ Isotype Control Recombinant Antibody

Catalog Number:CL594-98600



Basic Information

Catalog Number:

CL594-98600

Size:

100ug , 500 ug/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

252253A1

Recommended Dilutions:

FC: 0.25 ug per 10^6 cells in a 100 μ l suspension

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

n/a

Positive Controls:

FC : human PBMCs, rat splenocytes, THP-1 cells

Background Information

This product is an isotype control antibody, which is used to estimate the non-specific binding of target primary antibodies due to Fc receptor binding or other protein-protein interactions. An isotype control antibody should have the same immunoglobulin type and be used at the same concentration as the test antibody.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were unstained (left, blue in histogram) or surface stained with 0.25 ug CoraLite®594 Mouse IgG2a, λ Isotype Control RecAb (CL594-98600, Clone: 252253A1) (middle, red in histogram) or 0.25 ug CoraLite®594 Anti-Human CD8 (UCHT4) Mouse IgG2a Recombinant Antibody (CL594-65559, Clone: UCHT4) (right, green in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.

1×10^6 rat splenocytes were unstained (left, blue in histogram) or surface stained with 0.25 ug CoraLite®594 Mouse IgG2a, λ Isotype Control RecAb (CL594-98600, Clone: 252253A1) (middle, red in histogram) or 0.25 ug CoraLite®594 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (CL594-65607, Clone: W3/13) (right, green in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.