

For Research Use Only

CoraLite® Plus 647-conjugated MAGEA4 Polyclonal antibody



Catalog Number: CL647-12508

Basic Information

Catalog Number: CL647-12508	GenBank Accession Number: BC017723	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 4103	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Source: Rabbit	Full Name: melanoma antigen family A, 4	
Isotype: IgG	Calculated MW: 317 aa, 35 kDa	
Immunogen Catalog Number: AG3193	Observed MW: 45 kDa	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.2 ug CoraLite® Plus 647 Anti-Human MAGEA4 (CL647-12508) (red), or 0.2 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).