

For Research Use Only

# CoraLite® Plus 647 Anti-Mouse CSF1R/CD115 (AFS98)

Catalog Number: CL647-65095



## Basic Information

Catalog Number:

CL647-65095

Size:

100ug , 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

Excitation/Emission maxima  
wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> Balb/c mouse peritoneal macrophages were surface stained with 0.25 ug CoraLite® Plus 647 Anti-Mouse CD115 (c-fms) (AFS98) (CL647-65095, Clone: AFS98), or 0.25 ug CoraLite® Plus 647 Rat IgG2a Isotype Control (2A3) (CL647-65209, Clone: 2A3). Cells were not fixed.