

For Research Use Only

CoraLite® Plus 647 Anti-Human CD45 (HI30)



Catalog Number: **CL647-65109**

Basic Information

Catalog Number: CL647-65109	GenBank Accession Number: BC014239	Purification Method: Affinity purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 5788	CloneNo.: HI30
Source: Mouse	Full Name: protein tyrosine phosphatase, receptor type, C	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG1, kappa		Excitation/Emission maxima wavelengths: 654 nm / 674 nm

Applications

Tested Applications: FC, IF	Positive Controls: IF : Jurkat cells,
Species Specificity: Human	

Background Information

CD45, also known as protein tyrosine phosphatase, receptor type C, is a type I transmembrane protein expressed on the surface of all haematopoietic cells with the exception of erythrocytes and platelets (PMID: 3489673; 28615666). CD45 is a pan-haematopoietic cell marker and has been shown to be essential for T- and B-cell activation and signalling (PMID: 9429890; 16378097).


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 5 μ l CoraLite® Plus 647 Anti-Human CD45 (CL647-65109, Clone: HI30), or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using CoraLite® Plus 647 CD45 antibody (CL647-65109, Clone: HI30) at dilution of 1:200.