

For Research Use Only

CoraLite® Plus 647 Anti-Mouse TER-119 (TER-119)



Catalog Number: **CL647-65149**

Basic Information

Catalog Number:

CL647-65149

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

GeneID (NCBI):

104231

Full Name:

lymphocyte antigen 76

Purification Method:

Affinity purification

CloneNo.:

TER-119

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

TER-119 is a mouse erythroid lineage-specific monoclonal antibody that reacts with erythroid cells at differentiation stages from early proerythroblast to mature erythrocyte, but not with cells showing typical erythroid blast-forming unit (BFU-E) and erythroid colony-forming unit (CFU-E) activities (PMID: 1975515; 10848813). TER-119 recognizes a 52-kDa molecule on erythrocyte membranes (PMID: 10848813). TER-119 antigen is a molecule associated with cell-surface glycophorin A but not with glycophorin A itself (PMID: 10848813).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse bone marrow cells were surface co-stained with CoraLite® 488 Anti-Mouse CD45 (30-F11) and 0.25 ug CoraLite® Plus 647 Anti-Mouse TER-119 (CL647-65149, Clone:TER-119) or 0.25 ug CoraLite® Plus 647-conjugated Rat IgG2b isotype control. Cells were not fixed.