

For Research Use Only

CoraLite® Plus 647 Anti-Human CD1a (HI149)



Catalog Number: **CL647-65172**

Basic Information

Catalog Number:
CL647-65172

Size:
100tests , 5 µl/test

Source:
Mouse

Isotype:
IgG1, kappa

GenBank Accession Number:
BC031645

GeneID (NCBI):
909

Full Name:
CD1a molecule

Calculated MW:
327 aa, 37 kDa

Purification Method:
Purified by protein-A affinity chromatography

CloneNo.:
HI149

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ MOLT-4 cells were surface stained with 5 ul CoraLite® Plus 647 Anti-Human CD1a (CL647-65172, Clone:HI 149). Cells were not fixed.