

For Research Use Only

# CoraLite® Plus 647 Anti-Human STRO-1 (STRO-1)



Catalog Number: **CL647-65184**

## Basic Information

Catalog Number:  
**CL647-65184**

Size:  
100tests , 5 µl/test

Source:  
Mouse

Isotype:  
IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Purified by precipitation and chromatography

CloneNo.:

STRO-1

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human whole blood cells were surface stained with 5 ul CoraLite® Plus 647 Anti-Human STRO-1 (CL647-65184, Clone: STRO-1). Cells were not fixed. CD45dim, CD34+ cells were gated.