

For Research Use Only

CoraLite® Plus 647 Anti-Mouse CD41 (MWReg30)



Catalog Number: CL647-65189

Basic Information

Catalog Number: CL647-65189	GenBank Accession Number: BC120493	Purification Method: Affinity purification
Size: 100ug , 500 µg/ml	GeneID (NCBI): 16399	CloneNo.: MWReg30
Source: Rat	Full Name: integrin alpha 2b	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Isotype: IgG1, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Balb/c mouse peripheral blood platelets were surface stained with 0.5 ug CoraLite® Plus 647 Anti-Mouse CD41 (CL647-65189, Clone:MWRReg30), or 0.5 ug Isotype Control. Cells were not fixed.



1X10⁶ Balb/c mouse peripheral blood platelets were surface stained with 0.5 ug CoraLite® Plus 647 Anti-Mouse CD41 (CL647-65189, Clone:MWRReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.