

For Research Use Only

CoraLite® Plus 647 Anti-Mouse CD49b (HMA2)

Catalog Number: CL647-65245



Basic Information

Catalog Number: CL647-65245	GenBank Accession Number: BC065139	Purification Method: Affinity purification
Size: 100ug , 500 µg/ml	GeneID (NCBI): 16398	CloneNo.: HMA2
Source: Armenian Hamster	Full Name: integrin alpha 2	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Isotype: IgG		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse NK1.1 (CD161) (FITC-65138, Clone: PK136) and 0.5 ug CoraLite® Plus 647 Anti-Mouse CD49b (CL647-65245, Clone: HMA 2) or 0.5 ug CoraLite® Plus 647-conjugated Armenian Hamster IgG. Cells were not fixed.