

For Research Use Only

CoraLite® Plus 647 Anti-Human CD96 (NK92.39)



Catalog Number: CL647-65250

Basic Information

Catalog Number: CL647-65250	GenBank Accession Number: BC020749	Purification Method: Affinity purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 10225	CloneNo.: NK92.39
Source: Mouse	Full Name: CD96 molecule	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Isotype: IgG1, kappa	Calculated MW: 402 aa, 46 kDa	

Applications

Tested Applications:
FC
Species Specificity:
Human

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface co-stained with Anti-Human CD19 and 5 ul CoraLite® Plus 647 Anti-Human CD96 (CL647-65250, Clone: NK92.39), or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.