

For Research Use Only



# CoraLite® Plus 647 Anti-Rat CD43 (W3/13) Mouse IgG1 Recombinant Antibody

Catalog Number: **CL647-65535**

## Basic Information

Catalog Number:

CL647-65535

Size:

100ug , 500 µg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Affinity purification

CloneNo.:

W3/13

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

Rat

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 647 Anti-Rat CD43 (W3/13) Mouse Recombinant Antibody (CL647-65535, Clone:W3/13), or CoraLite® Plus 647 Mouse IgG1 Isotype Control (MOPC-21) (CL647-65124, Clone: MOPC-21). Cells were not fixed.