

For Research Use Only

CoraLite® Plus 647 Anti-Rat CD11b (OX42) Rabbit Recombinant Antibody



Catalog Number:CL647-65548

Basic Information

Catalog Number:

CL647-65548

Size:

100ug , 500 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_012711

GeneID (NCBI):

25021

Full Name:

integrin alpha M

Calculated MW:

127 kDa

Purification Method:

Affinity purification

CloneNo.:

OX42

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

Rat

Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11b, also known as Integrin alpha M or CR3A, belongs to the integrin alpha chain family. CD11b forms an α/β heterodimer with CD18 (integrin β 2). CD11b/CD18 is implicated in various adhesive interactions of monocytes, macrophages and granulocytes as well as in mediating the uptake of complement-coated particles and pathogens (PMID: 9558116; 20008295). CD11b/CD18 is a receptor for the complement protein fragment iC3b, and is also a receptor for fibrinogen, factor X and ICAM1 (PMID: 2971974; 15485828).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat bone marrow cells were intracellularly stained with 0.5 ug CoraLite® Plus 647 Anti-Rat CD11b (OX42) Rabbit Recombinant Antibody (CL647-65548, Clone:OX42), or CoraLite® Plus 647-conjugated Rabbit IgG control Rabbit PolyAb (CL647-30000). Cells were not fixed.