

For Research Use Only

CoraLite® Plus 647-conjugated TOM20 Recombinant monoclonal antibody

Catalog Number: CL647-80501



Basic Information

Catalog Number: CL647-80501	GenBank Accession Number: BC000882	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 9804	CloneNo.: 1D22
Source: Rabbit	UNIPROT ID: Q15388	Recommended Dilutions: IF/ICC: 1:50-1:500
Isotype: IgG	Full Name: translocase of outer mitochondrial membrane 20 homolog (yeast)	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Immunogen Catalog Number: AG2378	Calculated MW: 145 aa, 16 kDa	
	Observed MW: 16 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HepG2 cells,
Species Specificity: human, mouse, rat, pig, chicken, zebrafish	

Background Information

TOM20, also named as KIAA0016, belongs to the Tom20 family. It is a central component of the receptor complex responsible for the recognition and translocation of cytosolically synthesized mitochondrial preproteins. Together with TOM22, TOM20 functions as the transit peptide receptor at the surface of the mitochondrion outer membrane and facilitates the movement of preproteins into the TOM40 translocation pore. TOM20 is characterized as major docking receptors to mediate the recognition by different mechanisms.


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CoraLite® Plus 647 TOM20 antibody (CL647-80501, Clone: 1D22) at dilution of 1:200.