

For Research Use Only

CoraLite® Plus 647-conjugated transgelin/SM22 Recombinant monoclonal antibody

Catalog Number: CL647-83922-2



Basic Information

Catalog Number: CL647-83922-2	GenBank Accession Number: BC004927	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 6876	CloneNo.: 240864A11
Source: Rabbit	UNIPROT ID: Q01995	Recommended Dilutions: IF-P: 1:50-1:500
Isotype: IgG	Full Name: transgelin	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Immunogen Catalog Number: AG0764	Calculated MW: 22 kDa	
	Observed MW: 19-22 kDa	

Applications

Tested Applications: IF-P	Positive Controls: IF-P : mouse heart tissue,
Species Specificity: human, mouse, rat	

Background Information

The transgelin family is a protein group belonging to the 22kd actin-related calponin superfamily. Of all three isoforms, transgelin 1 is the best characterized. Transgelin 1, or SM22alpha, is a specific marker for differentiated smooth muscle cells. Transgelin 2, also known as SM22 beta, is expressed by both smooth and non-smooth muscle cells in a temporally and spatially regulated pattern. Transgelin 3, also known as NP25, is only found in highly differentiated neuronal cells. This antibody was generated against full-length transgelin 1 protein. It can cross-react with other two transgenes based on the sequence similarity.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse heart tissue using CoraLite® Plus 647 transgelin/SM22 antibody (CL647-83922-2, Clone: 240864A11) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).