

For Research Use Only

CoraLite® Plus 647-conjugated PCNP Recombinant monoclonal antibody

Catalog Number:CL647-86075



Basic Information

Catalog Number:

CL647-86075

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30951

GenBank Accession Number:

BC013916

GeneID (NCBI):

57092

UNIPROT ID:

Q8WW12

Full Name:

PEST proteolytic signal containing
nuclear protein

Calculated MW:

19 kDa

Observed MW:

24-28 kDa

Purification Method:

Protein A purification

CloneNo.:

250576B6

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HepG2 cells,

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using Coralite® Plus 647 PCNP antibody (CL647-86075, Clone: 250576B6) at dilution of 1:200, CL488-phalloidin (green).