

For Research Use Only

FcZero-rAb® CoraLite® Plus 647 Anti-Human FLT3/CD135 Rabbit Recombinant Antibody

Catalog Number: CL647-FcA98106



Basic Information

Catalog Number:

CL647-FcA98106

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC036028

GeneID (NCBI):

2322

UNIPROT ID:

P36888

Full Name:

fms-related tyrosine kinase 3

Calculated MW:

993 aa, 113 kDa

Purification Method:

Protein A purification

CloneNo.:

241414F8

Recommended Dilutions:

FC: 5 ul per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : Reh cells,

Background Information

FLT3 (also known as CD135 or FLK2) is a tyrosine-protein kinase that acts as a cell-surface receptor for the cytokine FLT3LG and regulates differentiation, proliferation, and survival of hematopoietic progenitor cells and of dendritic cells. FLT3 was originally identified by its expression in hematopoietic stem/progenitor cells (PMID: 7507245). It is important for the normal development of hematopoietic stem/progenitor cells. Mutations that result in the constitutive activation of this receptor result in acute myeloid leukemia and acute lymphoblastic leukemia.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ Reh cells (left) or U-937 cells (right) were surface stained with 5 μ l CoraLite[®] Plus 647 Anti-Human FLT3/CD135 Rabbit RecAb (CL647-FcA98106, Clone: 241414F8) (red) or FcZero-rAb[®] CoraLite[®] Plus 647 Rabbit IgG Isotype Control Recombinant Antibody (CL647-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.